US Probe Calibration with LSQRREcipes 1.0

Generated by Doxygen 1.8.0

Thu May 17 2012 14:33:36

Contents

1	Nam	nespace	Index									1
	1.1	Names	pace List					 	 	 	 	1
2	Clas	s Index										3
	2.1	Class	_ist					 	 	 	 	3
3	File	Index										5
	3.1	File Lis	t					 	 	 	 	5
4	Nam	nespace	Documentatio	n								7
	4.1	Ui Nan	nespace Refere	nce				 	 	 	 	7
5	Clas	s Docu	mentation									9
	5.1	Calibra	tion Class Refe	rence				 	 	 	 	9
		5.1.1	Detailed Desc	ription				 	 	 	 	10
		5.1.2	Member Typed	def Documen	itation .			 	 	 	 	10
			5.1.2.1 Data	аТуре				 	 	 	 	10
		5.1.3	Member Funct	ion Docume	ntation			 	 	 	 	10
			5.1.3.1 Cali	brate				 	 	 	 	10
			5.1.3.2 Clea	arlmagePoint	is			 	 	 	 	10
			5.1.3.3 Clea	arTransforma	itions .			 	 	 	 	10
			5.1.3.4 getE	EstimatedUS	Calibration	onParar	neters	 	 	 	 	10
			5.1.3.5 Inse	rtImagePoin	ts			 	 	 	 	10
			5.1.3.6 Inse	rtTransforma	ations .			 	 	 	 	10
			5.1.3.7 New					 	 	 	 	10
		5.1.4	Member Data	Documentati	ion			 	 	 	 	10
			5.1.4.1 data	ι				 	 	 	 	10
			5.1.4.2 estir	matedUSCali	ibrationF	Paramet	ers .	 	 	 	 	10
			5.1.4.3 imag	gePoints				 	 	 	 	10
			5.1.4.4 tran	sformations				 	 	 	 	11
	5.2	MainW	indow Class Re	ference				 	 	 	 	11
		5.2.1	Detailed Desc	ription				 	 	 	 	11
		5.2.2	Constructor &	Destructor D	ocumen	itation .		 	 	 	 	12

ii CONTENTS

	5.2.2.1	MainWindow	12
	5.2.2.2	~MainWindow	12
5.2.3	Member	Function Documentation	12
	5.2.3.1	addImages	12
	5.2.3.2	addLogText	12
	5.2.3.3	displaySelectedImage	12
	5.2.3.4	getDisplayWidget	12
	5.2.3.5	print	12
	5.2.3.6	probeCalibration	12
5.2.4	Member	Data Documentation	12
	5.2.4.1	Connections	12
	5.2.4.2	displayWidget	12
	5.2.4.3	imagesFilenames	12
	5.2.4.4	textOnTextArea	13
	5.2.4.5	ui	13
ProbeC	Calibration	Widget Class Reference	13
5.3.1	Detailed	Description	14
5.3.2	Construc	ctor & Destructor Documentation	14
	5.3.2.1	ProbeCalibrationWidget	14
	5.3.2.2	~ProbeCalibrationWidget	14
5.3.3	Member	Function Documentation	14
	5.3.3.1	calibrate	14
	5.3.3.2	crop	14
	5.3.3.3	cropProbeImage	14
	5.3.3.4	getCoordinates	15
	5.3.3.5	loadRotationsFile	15
	5.3.3.6	loadTranslationsFile	15
	5.3.3.7	saveCalibration	15
	5.3.3.8	setImage	15
	5.3.3.9	setImageStack	15
	5.3.3.10	setMainWindow	15
5.3.4	Member	Data Documentation	15
	5.3.4.1	calibrationParameters	15
	5.3.4.2	coords	15
	5.3.4.3	image	15
	5.3.4.4	imageStack	16
	5.3.4.5	mainWindow	16
	5.3.4.6	rotations	16
	5.3.4.7	rotations_2	16
	5.3.4.8	translations	16
	5.2.4 ProbeC 5.3.1 5.3.2	5.2.2.2 5.2.3 5.2.3.1 5.2.3.2 5.2.3.4 5.2.3.5 5.2.3.6 5.2.4.1 5.2.4.2 5.2.4.3 5.2.4.3 5.2.4.5 ProbeCalibration 5.3.1 Detailed 5.3.2 Construct 5.3.2.1 5.3.2.2 5.3.3 Member 5.3.3.1 5.3.3.2 5.3.3.3 5.3.3.4 5.3.3.5 5.3.3.6 5.3.3.6 5.3.3.7 5.3.3.8 5.3.3.9 5.3.3.10 5.3.4.1 5.3.4.2 5.3.4.1 5.3.4.2 5.3.4.3 5.3.4.2 5.3.4.3 5.3.4.4 5.3.4.5 5.3.4.6 5.3.4.5 5.3.4.6 5.3.4.7	5.2.3 Member Function Documentation 5.2.3.1 addImages 5.2.3.2 addLogText 5.2.3.3 displaySelectedImage 5.2.3.4 getDisplayWidget 5.2.3.5 probeCalibration 5.2.4 Member Data Documentation 5.2.4.1 Connections 5.2.4.2 displayWidget 5.2.4.3 imagesFilenames 5.2.4.4 taxOnTextArea 5.2.4.5 ul ProbeCalibrationWidget Class Reference 5.3.1 Detailed Description 5.3.2 ProbeCalibrationWidget 5.3.2 ~ProbeCalibrationWidget 5.3.3.1 calibrate 5.3.3.2 crop 5.3.3.3 cropProbeImage 5.3.3.4 getCoordinates 5.3.3.5 loadTranslationsFile 5.3.3.7 saveCalibration 5.3.3.8 setImage 5.3.3.1 calibrationPrameters 5.3.4.2 coords 5.3.4.3 image 5.3.4.4 imageStack 5.3.4.5 mainWindow

CONTENTS

		5.3.4.9	workWithStack	16				
5.4	QVTKI	mageWidg	get Class Reference	16				
	5.4.1	Detailed Description						
	5.4.2	Construc	tor & Destructor Documentation	18				
		5.4.2.1	QVTKImageWidget	18				
		5.4.2.2	~QVTKImageWidget	18				
	5.4.3	Member	Function Documentation	18				
		5.4.3.1	displayImage	18				
		5.4.3.2	displaySelectedImage	18				
		5.4.3.3	getImageDisplayedIndex	18				
		5.4.3.4	getImageHeigth	18				
		5.4.3.5	getImageSize	18				
		5.4.3.6	getImageStack	18				
		5.4.3.7	getImageType	19				
		5.4.3.8	getImageViewer	19				
		5.4.3.9	getImageWidth	19				
		5.4.3.10	getNumOfDimesions	19				
		5.4.3.11	getPixelType	19				
		5.4.3.12	getQVTKWidget	19				
		5.4.3.13	getXPicked	19				
		5.4.3.14	getYPicked	19				
		5.4.3.15	setAndDisplayImage	19				
		5.4.3.16	setAndDisplayImage	20				
		5.4.3.17	setAndDisplayMultipleImages	20				
		5.4.3.18	setAndDisplayMultipleImages	20				
		5.4.3.19	setImageProperties	20				
		5.4.3.20	setXPicked	20				
		5.4.3.21	setYPicked	20				
	5.4.4	Member	Data Documentation	20				
		5.4.4.1	cornerAnnotation	20				
		5.4.4.2	imageDisplayedIndex	21				
		5.4.4.3	imageHeight	21				
		5.4.4.4	imageStack	21				
		5.4.4.5	imageType	21				
		5.4.4.6	imageViewer	21				
		5.4.4.7	imageWidth	21				
		5.4.4.8	isImageStackLoaded	21				
		5.4.4.9	itklmage	21				
		5.4.4.10	numDimensions	21				
		5.4.4.11	pixelType	21				

iv CONTENTS

			5.4.4.12	qvtkWidget	 21
			5.4.4.13	enderer	 21
			5.4.4.14	enwin	 22
			5.4.4.15	gbltkImage	 22
			5.4.4.16	rtkImage	 22
			5.4.4.17	(Picked	 22
			5.4.4.18	Position	 22
			5.4.4.19	Picked	 22
			5.4.4.20	/Position	 22
	5.5	QVTKI	mageWidge	tCommand Class Reference	 22
		5.5.1	Detailed D	escription	 23
		5.5.2	Constructo	r & Destructor Documentation	 23
			5.5.2.1	QVTKImageWidgetCommand	 23
			5.5.2.2	~QVTKImageWidgetCommand	 23
		5.5.3	Member F	unction Documentation	 23
			5.5.3.1	Execute	 23
			5.5.3.2	New	 23
			5.5.3.3	SetAnnotation	 23
			5.5.3.4	SetImageWidget	 23
			5.5.3.5	SetPicker	 23
		5.5.4	Member D	ata Documentation	 23
			5.5.4.1	Annotation	 23
			5.5.4.2	mageWidget	 23
			5.5.4.3	Picker	 24
6	File	Documo	entation		25
	6.1			Reference	 25
	6.2			deference	25
	6.3			rence	25
		6.3.1		ocumentation	25
				main	25
	6.4	mainw		le Reference	26
	6.5			Reference	26
	6.6			idget.cpp File Reference	26
		6.6.1		ocumentation	26
			6.6.1.1	setCoordsSize	 26
	6.7	Probe(CalibrationW	idget.h File Reference	 27
	6.8	QVTKI	mageWidge	t.cpp File Reference	 27
	6.9	QVTKI	mageWidge	t.h File Reference	 28
		6.9.1	Typedef Do	ocumentation	 28

CONTENTS		7

6.9.1.1	ImageType	28
6.9.1.2	RGBImageType	28
6.9.1.3	RGBPixelType	28
6.10 QVTKImageWide	getCommand.cpp File Reference	29
6.11 QVTKImageWide	getCommand.h File Reference	29

Chapter 1

Namespace Index

1.1	Namespace List
Here	is a list of all namespaces with brief descriptions:
1.6	

2 Namespace Index

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

Calibration
Implements LSQRRecepies methods
MainWindow
Main Window for user Interaction
ProbeCalibrationWidget
Obtain data for calibration process
QVTKImageWidget
Display VTK images
QVTKImageWidgetCommand
Interaction with mouse

Class Index

Chapter 3

File Index

3.1 File List

Here is a list of all files with brief descriptions:

Calibration.cpp	25
Calibration.h	25
main.cpp	25
mainwindow.cpp	26
mainwindow.h	26
ProbeCalibrationWidget.cpp	26
ProbeCalibrationWidget.h	27
QVTKImageWidget.cpp	27
QVTKImageWidget.h	28
QVTKImageWidgetCommand.cpp	29
QVTKImageWidgetCommand.h	29

6 File Index

Chapter 4

Namespace Documentation

4.1 Ui Namespace Reference

Names	pace	Do	cu	me	nta	tic	n

Chapter 5

Class Documentation

5.1 Calibration Class Reference

Implements LSQRRecepies methods.

```
#include <Calibration.h>
```

Public Types

typedef
 lsqrRecipes::SingleUnknownPointTargetUSCalibrationParametersEstimator::DataType DataType

Public Member Functions

- void InsertTransformations (vnl_matrix< double > rotationMatrix, vnl_vector< double > translation) insert the rotation matrix of an image to transformations
- void ClearTransformations ()

clear transformations

void InsertImagePoints (double p[2])

insert the crosswire point of an image to imagePoints

void ClearImagePoints ()

clear imagePoints

• bool Calibrate ()

estimate calibration parameters

 $\bullet \ \ std:: vector < double > getEstimatedUSCalibrationParameters \ () \\$

Static Public Member Functions

static Calibration * New ()
 Constructor of the class.

Private Attributes

- std::vector< lsqrRecipes::Frame > transformations
- std::vector< lsqrRecipes::Point2D > imagePoints

contains the crosswire point in all images

std::vector< DataType > data

contain the data of all images

• std::vector< double > estimatedUSCalibrationParameters

5.1.1 Detailed Description

Implements LSQRRecepies methods.

contains the crosswire point in all images

This classs have the calibration methods implemented in LSQRRecipes to calibrate an Ultra Sound Probe with a cross wire phantom.

```
Member Typedef Documentation
5.1.2
       type def\ Is qr Recipes:: Single Unknown Point Target USC alibration Parameters Estimator:: Data Type
       Calibration::DataType
5.1.3 Member Function Documentation
5.1.3.1 bool Calibration::Calibrate ( )
estimate calibration parameters
5.1.3.2 void Calibration::ClearImagePoints ( )
clear imagePoints
5.1.3.3 void Calibration::ClearTransformations()
clear transformations
5.1.3.4 std::vector< double > Calibration::getEstimatedUSCalibrationParameters ( )
5.1.3.5 void Calibration::InsertImagePoints (double p[2])
insert the crosswire point of an image to imagePoints
5.1.3.6 void Calibration::InsertTransformations ( vnl_matrix< double > rotationMatrix, vnl_vector< double >
       translation )
insert the rotation matrix of an image to transformations
5.1.3.7 static Calibration* Calibration::New() [inline, static]
Constructor of the class.
5.1.4
       Member Data Documentation
5.1.4.1 std::vector < DataType > Calibration::data [private]
contain the data of all images
5.1.4.2 std::vector<double> Calibration::estimatedUSCalibrationParameters [private]
5.1.4.3 std::vector<lsqrRecipes::Point2D> Calibration::imagePoints [private]
```

5.1.4.4 std::vector<lsqrRecipes::Frame> Calibration::transformations [private]

The documentation for this class was generated from the following files:

- · Calibration.h
- · Calibration.cpp

5.2 MainWindow Class Reference

Main Window for user Interaction.

```
#include <mainwindow.h>
```

Public Member Functions

- MainWindow (QWidget *parent=0)
- ∼MainWindow ()
- void addLogText (QString str)

Writes at the logger.

QVTKImageWidget * getDisplayWidget ()

return this display widget

Private Slots

· void addImages ()

Add image folder to application.

• void displaySelectedImage (int idx)

Display selected image with the image slider.

void probeCalibration ()

Implements a ultrasound 3D probe calibration, for navigate with the probe.

• void print ()

Print the text on logger.

Private Attributes

• Ui::MainWindow * ui

main window object

QStringList imagesFilenames

The filename of each selected image.

QString textOnTextArea

Text on the logger.

- QVTKImageWidget * displayWidget
- vtkSmartPointer
 - < vtkEventQtSlotConnect > Connections

Connects with the probe calibration widget.

5.2.1 Detailed Description

Main Window for user Interaction.

This class is the main window of ht US Probe Calibration. Here the user can open the calibration images and open the calibration widget

5.2.2 Constructor & Destructor Documentation

5.2.2.1 MainWindow::MainWindow(QWidget* parent = 0) [explicit]

5.2.2.2 MainWindow::~MainWindow()

5.2.3 Member Function Documentation

5.2.3.1 void MainWindow::addImages() [private, slot]

Add image folder to application.

5.2.3.2 void MainWindow::addLogText (QString str)

Writes at the logger.

5.2.3.3 void MainWindow::displaySelectedImage (int *idx*) [private, slot]

Display selected image with the image slider.

5.2.3.4 QVTKImageWidget * MainWindow::getDisplayWidget ()

return this display widget

Parameters

out	this	display widget

5.2.3.5 void MainWindow::print() [private, slot]

Print the text on logger.

5.2.3.6 void MainWindow::probeCalibration() [private, slot]

Implements a ultrasound 3D probe calibration, for navigate with the probe.

5.2.4 Member Data Documentation

5.2.4.1 vtkSmartPointer<vtkEventQtSlotConnect> MainWindow::Connections [private]

Connects with the probe calibration widget.

5.2.4.2 QVTKImageWidget* **MainWindow::displayWidget** [private]

Central widget for display image purposes

5.2.4.3 QStringList MainWindow::imagesFilenames [private]

The filename of each selected image.

5.2.4.4 QString MainWindow::textOnTextArea [private]

Text on the logger.

5.2.4.5 Ui::MainWindow* MainWindow::ui [private]

main window object

The documentation for this class was generated from the following files:

- · mainwindow.h
- · mainwindow.cpp

5.3 ProbeCalibrationWidget Class Reference

Obtain data for calibration process.

```
#include <ProbeCalibrationWidget.h>
```

Public Member Functions

- ProbeCalibrationWidget (QWidget *parent=0)
- virtual ∼ProbeCalibrationWidget ()
- void setImageStack (std::vector< vtkSmartPointer< vtkImageData > > imageStack)

Set this stack of vtklmageData.

void setImage (vtkSmartPointer< vtkImageData > image)

Set this vtklmageData.

void setMainWindow (MainWindow *mainwindow)

Private Slots

• void crop ()

crop the images to delete the extra information

• void getCoordinates ()

get the cross point coordinates

· void calibrate ()

Calls Calibrate.h to estimate the calibration parameters.

void loadRotationsFile ()

Load the rotations file.

· void loadTranslationsFile ()

Load the translation file.

• void saveCalibration ()

Save the Estimated Parameters in a .txt file.

Private Member Functions

 vtkSmartPointer< vtkImageData > cropProbeImage (vtkSmartPointer< vtkImageData > image, int depth-Type)

Private Attributes

```
· bool workWithStack
```

• std::vector< vtkSmartPointer

< vtkImageData > > imageStack

an Array of vtklmageData to work

vnl matrix< double > translations

a vnl_matrix to store the translations of each image given by the tracker

vnl_matrix< double > rotations

a vnl_matrix to store the rotations of each image given by the tracker

- float rotations_2 [7][4]
- MainWindow * mainWindow
- vtkSmartPointer< vtkImageData > image

the vtkImageData to work

vnl_matrix< double > coords

a vnl_matrix to store the selected coordinates on each image

std::vector< double > calibrationParameters

the estimate calibration parameters by Calibration.h

5.3.1 Detailed Description

Obtain data for calibration process.

This class obtain the necesary data to use the Calibrtion class. It let the user to set the croos wire point on the images, load Rotation and Translation Data. It also allows the user to save the calibration estimated paraeters in a .txt file

5.3.2 Constructor & Destructor Documentation

5.3.2.1 ProbeCalibrationWidget::ProbeCalibrationWidget (QWidget * parent = 0)

Constructor

- **5.3.2.2 ProbeCalibrationWidget::**~ProbeCalibrationWidget() [virtual]
- 5.3.3 Member Function Documentation
- **5.3.3.1 void ProbeCalibrationWidget::calibrate()** [private, slot]

Calls Calibrate.h to estimate the calibration parameters.

5.3.3.2 void ProbeCalibrationWidget::crop() [private, slot]

crop the images to delete the extra information

5.3.3.3 vtkSmartPointer< vtkImageData > ProbeCalibrationWidget::cropProbeImage (vtkSmartPointer< vtkImageData > image, int depthType) [private]

Crop ultrasound image depnding of the depth type

5.3.3.4 void ProbeCalibrationWidget::getCoordinates() [private, slot] get the cross point coordinates **5.3.3.5 void ProbeCalibrationWidget::loadRotationsFile()** [private, slot] Load the rotations file. **5.3.3.6 void ProbeCalibrationWidget::loadTranslationsFile()** [private, slot] Load the translation file. **5.3.3.7 void ProbeCalibrationWidget::saveCalibration()** [private, slot] Save the Estimated Parameters in a .txt file. 5.3.3.8 void ProbeCalibrationWidget::setImage (vtkSmartPointer< vtkImageData > image) Set this vtkImageData. **Parameters** a smart Pointer of vtklmageData in 5.3.3.9 void ProbeCalibrationWidget::setImageStack (std::vector< vtkSmartPointer< vtkImageData >> imageStack Set this stack of vtkImageData. **Parameters** in a std Vector of vtkImageData 5.3.3.10 void ProbeCalibrationWidget::setMainWindow (MainWindow * mainwindow) Set the window to display the crop images 5.3.4 Member Data Documentation **5.3.4.1** std::vector<double> ProbeCalibrationWidget::calibrationParameters [private] the estimate calibration parameters by Calibration.h **5.3.4.2** vnl_matrix<double> ProbeCalibrationWidget::coords [private] a vnl matrix to store the selected coordinates on each image

the vtkImageData to work

5.3.4.3 vtkSmartPointer<vtkImageData> ProbeCalibrationWidget::image [private]

```
5.3.4.4 std::vector< vtkSmartPointer<vtkImageData>> ProbeCalibrationWidget::imageStack [private]
an Array of vtkImageData to work

5.3.4.5 MainWindow* ProbeCalibrationWidget::mainWindow [private]
the main window to call it

5.3.4.6 vnl_matrix<double> ProbeCalibrationWidget::rotations [private]
a vnl_matrix to store the rotations of each image given by the tracker

5.3.4.7 float ProbeCalibrationWidget::rotations_2[7][4] [private]

5.3.4.8 vnl_matrix<double> ProbeCalibrationWidget::translations [private]
a vnl_matrix to store the translations of each image given by the tracker

5.3.4.9 bool ProbeCalibrationWidget::workWithStack [private]
```

The documentation for this class was generated from the following files:

- · ProbeCalibrationWidget.h
- ProbeCalibrationWidget.cpp

5.4 QVTKImageWidget Class Reference

```
Display VTK images.
```

#include <QVTKImageWidget.h>

Public Member Functions

- QVTKImageWidget (QWidget *parent=0)
- virtual ~QVTKImageWidget ()
- void setAndDisplayImage (QString imageFilename)
- $\bullet \ \ void \ set And Display Image \ (vtk Smart Pointer < vtk Image Data > image) \\$

Sets and display the given vtklmageData.

void setAndDisplayMultipleImages (QStringList filenames)

Set and display multiple images from a given images filenames. Display the image corresponding to the first element on the filenmaes list.

void setAndDisplayMultipleImages (std::vector< vtkSmartPointer< vtkImageData > > imageStack)

Set and display multiple images from a given vtklmageData Array.

void displaySelectedImage (int idx)

display an image stored in this imageStack.

• QString getPixelType ()

Returns the pixel type in loaded images.

• QString getImageType ()

Returns the type of image displayed.

• QString getNumOfDimesions ()

Return the numer of dimensions of the image.

- std::vector< vtkSmartPointer
 - < vtkImageData > > getImageStack ()

return this image stack

- int * getImageSize ()
- int getImageWidth ()
- · int getImageHeigth ()
- int getXPicked ()

Return the mouse x coordinate position when mouse left button is pressed.

• int getYPicked ()

Return the mouse y coordinate position when mouse left button is pressed.

void setXPicked (int xPosition)

Set the mouse x coordinate position when mouse left button is pressed.

void setYPicked (int yPosition)

Set the mouse y coordinate position when mouse left button is pressed.

vtkSmartPointer< vtkImageViewer2 > getImageViewer ()

Return this widget image viewer.

QVTKWidget * getQVTKWidget ()

Return this qvtkWidget.

• int getImageDisplayedIndex ()

If an image stack is loaded, then return the index in the image stack of displayed image.

Public Attributes

• bool isImageStackLoaded

Flag to know if it's displayed an image stack.

Private Member Functions

- void setImageProperties (bool verbose)
- void displayImage (vtkImageData *image)

Private Attributes

- QVTKWidget * qvtkWidget
- ImageType::Pointer itkImage
- RGBImageType::Pointer rgbltkImage
- vtkSmartPointer< vtkImageData > vtkImage
- vtkSmartPointer< vtkRenderWindow > renwin
- vtkSmartPointer< vtkRenderer > renderer
- std::vector< vtkSmartPointer
 - < vtkImageData > > imageStack

A vtkImageData Vector for keep the image references when load an image stack.

- std::string pixelType
- int imageType
- size_t numDimensions
- int imageWidth
- · int imageHeight
- int xPosition
- int yPosition
- int xPicked
- · int yPicked
- int imageDisplayedIndex

- vtkSmartPointer< vtkImageViewer2 > imageViewer
- vtkSmartPointer
 - < vtkCornerAnnotation > cornerAnnotation

Object for display information in the corners of the vtkImageViewer2.

5.4.1 Detailed Description

Display VTK images.

This class allows the user to display vtkUsers on QT

5.4.2 Constructor & Destructor Documentation

5.4.2.1 QVTKImageWidget::QVTKImageWidget (QWidget * parent = 0)

Constructor for this ImageWidget

5.4.2.2 QVTKImageWidget::~QVTKImageWidget() [virtual]

Destructor

5.4.3 Member Function Documentation

5.4.3.1 void QVTKImageWidget::displayImage(vtkImageData * image) [private]

Display the given vtkImage

5.4.3.2 void QVTKImageWidget::displaySelectedImage (int idx)

display an image stored in this imageStack.

Parameters

in	the	index in the stack position of the image

5.4.3.3 int QVTKImageWidget::getImageDisplayedIndex ()

If an image stack is loaded, then return the index in the image stack of displayed image.

5.4.3.4 int QVTKImageWidget::getImageHeigth()

returns this image heigth

5.4.3.5 int* QVTKImageWidget::getImageSize()

returns an array with the width and height of the image

 $5.4.3.6 \quad std:: vector < vtkSmartPointer < vtkImageData > > QVTKImageWidget:: getImageStack (\quad)$

return this image stack

5.4.3.7 QString QVTKImageWidget::getImageType()

Returns the type of image displayed.

5.4.3.8 vtkSmartPointer< vtkImageViewer2 > QVTKImageWidget::getImageViewer()

Return this widget image viewer.

Parameters

out	imageViewer	vtkImageViewer2 target 2D image.

5.4.3.9 int QVTKImageWidget::getImageWidth()

returns this image width

5.4.3.10 QString QVTKImageWidget::getNumOfDimesions()

Return the numer of dimensions of the image.

5.4.3.11 QString QVTKImageWidget::getPixelType()

Returns the pixel type in loaded images.

5.4.3.12 QVTKWidget * QVTKImageWidget::getQVTKWidget ()

Return this qvtkWidget.

Parameters

out	the	QVTKWidget

5.4.3.13 int QVTKImageWidget::getXPicked ()

Return the mouse x coordinate position when mouse left button is pressed.

Parameters

out	int	x position
-----	-----	------------

5.4.3.14 int QVTKImageWidget::getYPicked ()

Return the mouse y coordinate position when mouse left button is pressed.

Parameters

out	int	y position

5.4.3.15 void QVTKImageWidget::setAndDisplayImage (QString imageFilename)

Sets and display an image from a given image path

5.4.3.16 void QVTKImageWidget::setAndDisplayImage (vtkSmartPointer< vtkImageData > image)

Sets and display the given vtkImageData.

Parameters

lin	l a	vtklmageData to set and display
	_ ~	Thim age 2 and to cot and anopia,

5.4.3.17 void QVTKImageWidget::setAndDisplayMultipleImages (QStringList filenames)

Set and display multiple images from a given images filenames. Display the image corresponding to the first element on the filenmaes list.

Parameters

in	а	QStringList that contain the filename of each image
----	---	---

5.4.3.18 void QVTKImageWidget::setAndDisplayMultipleImages (std::vector< vtkSmartPointer< vtkImageData > imageStack)

Set and display multiple images from a given vtkImageData Array.

Parameters

in	а	std::vector of vtkImageData
	_ ~	otanioto o rimmagozata

5.4.3.19 void QVTKImageWidget::setImageProperties (bool *verbose* **)** [private]

Set the needed image properties (pixelType, imageType, num of dimensions)

5.4.3.20 void QVTKImageWidget::setXPicked (int xPosition)

Set the mouse x coordinate position when mouse left button is pressed.

Parameters

out	int	x position

5.4.3.21 void QVTKImageWidget::setYPicked (int yPosition)

Set the mouse y coordinate position when mouse left button is pressed.

Parameters

out	int	y position

5.4.4 Member Data Documentation

5.4.4.1 vtkSmartPointer<vtkCornerAnnotation> QVTKImageWidget::cornerAnnotation [private]

Object for display information in the corners of the vtklmageViewer2.

```
5.4.4.2 int QVTKImageWidget::imageDisplayedIndex [private]
If image stack is displayed this sets a reference to current image displayed
5.4.4.3 int QVTKImageWidget::imageHeight [private]
Heigth of the image
5.4.4.4 std::vector< vtkSmartPointer<vtkImageData> > QVTKImageWidget::imageStack [private]
A vtkImageData Vector for keep the image references when load an image stack.
5.4.4.5 int QVTKImageWidget::imageType [private]
the number of scalar components in the image 1 => grayscale, 3 => rgb
5.4.4.6 vtkSmartPointer<vtkImageViewer2> QVTKImageWidget::imageViewer [private]
the image viewer for display images
5.4.4.7 int QVTKImageWidget::imageWidth [private]
Width of the image
5.4.4.8 bool QVTKImageWidget::isImageStackLoaded
Flag to know if it's displayed an image stack.
5.4.4.9 ImageType::Pointer QVTKImageWidget::itkImage [private]
The grayscale image displayed in this widget
5.4.4.10 size_t QVTKImageWidget::numDimensions [private]
The number of the image dimensions
5.4.4.11 std::string QVTKImageWidget::pixelType [private]
The type of the image pixels
5.4.4.12 QVTKWidget* QVTKImageWidget::qvtkWidget [private]
The QVTKWidget for display and interact with the images
5.4.4.13 vtkSmartPointer<vtkRenderer> QVTKImageWidget::renderer [private]
The VTK renderer
```

- QVTKImageWidget.h
- QVTKImageWidget.cpp

5.5 QVTKImageWidgetCommand Class Reference

Interaction with mouse.

#include <QVTKImageWidgetCommand.h>

Public Member Functions

- QVTKImageWidgetCommand ()
- ∼QVTKImageWidgetCommand ()
- void SetPicker (vtkSmartPointer< vtkPropPicker > picker)
- void SetAnnotation (vtkSmartPointer< vtkCornerAnnotation > annotation)
- void SetImageWidget (QVTKImageWidget *imageWidget)

Set the 2d image widget related to this 2d event manager.

virtual void Execute (vtkObject *, unsigned long vtkNotUsed(event), void *)

Static Public Member Functions

static QVTKImageWidgetCommand * New ()

Private Attributes

- vtkSmartPointer< vtkPropPicker > Picker
- vtkSmartPointer
 - < vtkCornerAnnotation > Annotation
- QVTKImageWidget * ImageWidget

5.5.1 Detailed Description

Interaction with mouse.

The mouse motion callback, to pick the image and recover pixel values

- 5.5.2 Constructor & Destructor Documentation
- 5.5.2.1 QVTKImageWidgetCommand::QVTKImageWidgetCommand()
- 5.5.2.2 QVTKImageWidgetCommand::~QVTKImageWidgetCommand()
- 5.5.3 Member Function Documentation
- 5.5.3.1 void QVTKImageWidgetCommand::Execute (vtkObject * , unsigned long vtkNotUsedevent, void *) [virtual]
- 5.5.3.2 QVTKImageWidgetCommand * QVTKImageWidgetCommand::New() [static]
- 5.5.3.3 void QVTKImageWidgetCommand::SetAnnotation (vtkSmartPointer< vtkCornerAnnotation > annotation)
- 5.5.3.4 void QVTKImageWidgetCommand::SetImageWidget (QVTKImageWidget * imageWidget)

Set the 2d image widget related to this 2d event manager.

Parameters

in	viewer	QVTKImageWidget target 2D image

- 5.5.3.5 void QVTKImageWidgetCommand::SetPicker (vtkSmartPointer< vtkPropPicker > picker)
- 5.5.4 Member Data Documentation
- 5.5.4.1 vtkSmartPointer<vtkCornerAnnotation> QVTKImageWidgetCommand::Annotation [private]

Pointer to the annotation

5.5.4.2 QVTKImageWidget* QVTKImageWidgetCommand::ImageWidget [private]

The widget related to the mouse events

5.5.4.3 vtkSmartPointer<vtkPropPicker> QVTKImageWidgetCommand::Picker [private]

Pointer to the picker

The documentation for this class was generated from the following files:

- QVTKImageWidgetCommand.h
- QVTKImageWidgetCommand.cpp

Chapter 6

File Documentation

6.1 Calibration.cpp File Reference

```
#include "Calibration.h"
```

6.2 Calibration.h File Reference

```
#include "SinglePointTargetUSCalibrationParametersEstimator.h"
#include "RANSAC.h"
```

Classes

• class Calibration

Implements LSQRRecepies methods.

6.3 main.cpp File Reference

```
#include <QtGui/QApplication>
#include "mainwindow.h"
```

Functions

int main (int argc, char *argv[])

Main function. Creates a new MainWindow object.

6.3.1 Function Documentation

```
6.3.1.1 int main ( int argc, char * argv[] )
```

Main function. Creates a new MainWindow object.

26 File Documentation

6.4 mainwindow.cpp File Reference

```
#include "ui_mainwindow.h"
#include "mainwindow.h"
#include "ProbeCalibrationWidget.h"
#include <QVBoxLayout>
#include <vtkEventQtSlotConnect.h>
```

6.5 mainwindow.h File Reference

```
#include <QMainWindow>
#include <QtGui>
#include "QVTKImageWidget.h"
#include <vtkEventQtSlotConnect.h>
```

Classes

· class MainWindow

Main Window for user Interaction.

Namespaces

• namespace Ui

6.6 ProbeCalibrationWidget.cpp File Reference

```
#include "ProbeCalibrationWidget.h"
#include "Calibration.h"
#include <QErrorMessage>
#include <QString>
#include <QFile>
#include <QTextStream>
#include <vtkExtractVOI.h>
#include <vnl/vnl_quaternion.h>
#include <vnl/vnl_vector_fixed.h>
#include <vnl/algo/vnl_levenberg_marquardt.h>
#include <vnl/vnl_double_2.h>
```

Variables

• bool setCoordsSize = false

6.6.1 Variable Documentation

6.6.1.1 bool setCoordsSize = false

6.7 ProbeCalibrationWidget.h File Reference

```
#include "ui_ProbeCalibrationWidget.h"
#include "mainwindow.h"
#include <QWidget>
#include <vtkSmartPointer.h>
#include <vtkImageData.h>
#include <vnl/vnl_matrix.h>
#include <string>
#include <fstream>
#include <stdio.h>
```

Classes

• class ProbeCalibrationWidget

Obtain data for calibration process.

6.8 QVTKImageWidget.cpp File Reference

```
#include "QVTKImageWidget.h"
#include "QVTKImageWidgetCommand.h"
#include <QSize.h>
#include <QBoxLayout>
#include <QString>
#include <itkImage.h>
#include <itkImageFileReader.h>
#include <vtkImageReader2.h>
#include <vtkImageReader2Factory.h>
#include <vtkCornerAnnotation.h>
#include <vtkTextProperty.h>
#include <vtkImageActor.h>
#include <vtkImageActor.h>
#include <vtkImageActor.h>
#include <vtkImageActor.h>
#include <vtkInteractorStyleImage.h>
```

28 File Documentation

6.9 QVTKImageWidget.h File Reference

```
#include <QtGui>
#include <QWidget>
#include <QVTKWidget.h>
#include <itkImage.h>
#include <itkRGBPixel.h>
#include <vtkSmartPointer.h>
#include <vtkImageData.h>
#include <vtkRenderWindow.h>
#include <vtkRenderer.h>
#include <vtkCamera.h>
#include <vtkImageActor.h>
#include <vtkCommand.h>
#include <vtkImageViewer2.h>
#include <vtkCornerAnnotation.h>
#include <vtkTransform.h>
#include <vnl/vnl_quaternion.h>
#include <vnl/vnl_matrix.h>
#include <vnl/vnl vector.h>
```

Classes

· class QVTKImageWidget

Display VTK images.

Typedefs

- typedef itk::RGBPixel
 unsigned char > RGBPixelType
- typedef itk::Image< unsigned char > ImageType
- typedef itk::Image
 RGBPixelType, 2 > RGBImageType

6.9.1 Typedef Documentation

- $\textbf{6.9.1.2} \quad \textbf{typedef itk::} \\ \textbf{lmage} \\ < \\ \textbf{RGBPixelType}, \\ \textbf{2} \\ > \\ \textbf{RGBImageType}$
- 6.9.1.3 typedef itk::RGBPixel < unsigned char > RGBPixelType

6.10 QVTKImageWidgetCommand.cpp File Reference

```
#include "QVTKImageWidgetCommand.h"
#include <vtkImageActor.h>
#include <vtkImageData.h>
#include <vtkInteractorStyleImage.h>
#include <vtkRenderWindow.h>
#include <vtkRenderWindowInteractor.h>
#include <vtkVariant.h>
#include <vtkMath.h>
#include <vtkCommand.h>
#include <vtkImageViewer2.h>
```

6.11 QVTKImageWidgetCommand.h File Reference

```
#include "QVTKImageWidget.h"
#include <vtkCommand.h>
#include <vtkPropPicker.h>
#include <vtkCornerAnnotation.h>
#include <vtkSmartPointer.h>
```

Classes

class QVTKImageWidgetCommand

Interaction with mouse.

Index

\sim MainWindow	ProbeCalibrationWidget, 14
MainWindow, 12	
\sim ProbeCalibrationWidget	data
ProbeCalibrationWidget, 14	Calibration, 10
~QVTKImageWidget	DataType
QVTKImageWidget, 18	Calibration, 10
~QVTKImageWidgetCommand	displayImage
QVTKImageWidgetCommand, 23	QVTKImageWidget, 18
a · · · · · · · · · · · · · · · · · · ·	displaySelectedImage
addImages	MainWindow, 12
MainWindow, 12	QVTKImageWidget, 18
addLogText	displayWidget
MainWindow, 12	MainWindow, 12
Annotation	
QVTKImageWidgetCommand, 23	estimatedUSCalibrationParameters
QV Milliage Widge to ominand, 20	Calibration, 10
Calibrate	Execute
Calibration, 10	QVTKImageWidgetCommand, 23
calibrate	
ProbeCalibrationWidget, 14	getCoordinates
Calibration, 9	ProbeCalibrationWidget, 14
	getDisplayWidget
Calibrate, 10	MainWindow, 12
ClearImagePoints, 10	getEstimatedUSCalibrationParameters
ClearTransformations, 10	Calibration, 10
data, 10	getImageDisplayedIndex
DataType, 10	QVTKImageWidget, 18
estimatedUSCalibrationParameters, 10	getImageHeigth
getEstimatedUSCalibrationParameters, 10	QVTKImageWidget, 18
imagePoints, 10	getImageSize
InsertImagePoints, 10	QVTKImageWidget, 18
InsertTransformations, 10	getImageStack
New, 10	QVTKImageWidget, 18
transformations, 10	getImageType
Calibration.cpp, 25	QVTKImageWidget, 18
Calibration.h, 25	getImageViewer
calibrationParameters	QVTKImageWidget, 19
ProbeCalibrationWidget, 15	getImageWidth
ClearImagePoints	QVTKImageWidget, 19
Calibration, 10	getNumOfDimesions
ClearTransformations	QVTKImageWidget, 19
Calibration, 10	getPixelType
Connections	QVTKImageWidget, 19
MainWindow, 12	
coords	getQVTKWidget
ProbeCalibrationWidget, 15	QVTKImageWidget, 19
cornerAnnotation	getXPicked
QVTKImageWidget, 20	QVTKImageWidget, 19
crop	getYPicked
ProbeCalibrationWidget, 14	QVTKImageWidget, 19
cropProbeImage	imaga
oropi robellilage	image

INDEX 31

ProbeCalibrationWidget, 15	New
imageDisplayedIndex	Calibration, 10
QVTKImageWidget, 20	QVTKImageWidgetCommand, 23
imageHeight	numDimensions
QVTKImageWidget, 21	QVTKImageWidget, 21
imagePoints	Picker
Calibration, 10	
imageStack	QVTKImageWidgetCommand, 23
ProbeCalibrationWidget, 15	pixelType
QVTKImageWidget, 21	QVTKImageWidget, 21
ImageType	print
QVTKImageWidget.h, 28	MainWindow, 12
imageType	probeCalibration
QVTKImageWidget, 21	MainWindow, 12
imageViewer	ProbeCalibrationWidget, 13
QVTKImageWidget, 21	\sim ProbeCalibrationWidget, 14
ImageWidget	calibrate, 14
QVTKImageWidgetCommand, 23	calibrationParameters, 15
imageWidth	coords, 15
QVTKImageWidget, 21	crop, 14
imagesFilenames	cropProbelmage, 14
MainWindow, 12	getCoordinates, 14
	image, 15
InsertImagePoints	imageStack, 15
Calibration, 10	loadRotationsFile, 15
InsertTransformations	
Calibration, 10	loadTranslationsFile, 15
isImageStackLoaded	mainWindow, 16
QVTKImageWidget, 21	ProbeCalibrationWidget, 14
itkImage	ProbeCalibrationWidget, 14
QVTKImageWidget, 21	rotations, 16
	rotations_2, 16
loadRotationsFile	saveCalibration, 15
ProbeCalibrationWidget, 15	setImage, 15
loadTranslationsFile	setImageStack, 15
ProbeCalibrationWidget, 15	setMainWindow, 15
	translations, 16
main	workWithStack, 16
main.cpp, 25	ProbeCalibrationWidget.cpp, 26
main.cpp, 25	setCoordsSize, 26
main, 25	ProbeCalibrationWidget.h, 27
MainWindow, 11	,
\sim MainWindow, 12	QVTKImageWidget, 16
addImages, 12	~QVTKImageWidget, 18
addLogText, 12	cornerAnnotation, 20
Connections, 12	displayImage, 18
displaySelectedImage, 12	displaySelectedImage, 18
displayWidget, 12	getImageDisplayedIndex, 18
getDisplayWidget, 12	getImageHeigth, 18
imagesFilenames, 12	getImageSize, 18
MainWindow, 12	getImageStack, 18
MainWindow, 12	getImageType, 18
print, 12	getImageViewer, 19
probeCalibration, 12	getImageWidth, 19
textOnTextArea, 12	getNumOfDimesions, 19
ui, 13	getPixelType, 19
mainWindow	getQVTKWidget, 19
ProbeCalibrationWidget, 16	getXPicked, 19
mainwindow.cpp, 26	getYPicked, 19
mainwindow.h, 26	imageDisplayedIndex, 20

32 INDEX

imageHeight, 21	rotations_2
imageStack, 21	ProbeCalibrationWidget, 16
imageType, 21	
imageViewer, 21	saveCalibration
imageWidth, 21	ProbeCalibrationWidget, 15
•	setAndDisplayImage
isImageStackLoaded, 21	QVTKImageWidget, 19
itklmage, 21	
numDimensions, 21	setAndDisplayMultipleImages
pixelType, 21	QVTKImageWidget, 20
QVTKImageWidget, 18	SetAnnotation
qvtkWidget, 21	QVTKImageWidgetCommand, 23
QVTKImageWidget, 18	setCoordsSize
renderer, 21	ProbeCalibrationWidget.cpp, 26
renwin, 21	setImage
rgbltkImage, 22	ProbeCalibrationWidget, 15
setAndDisplayImage, 19	setImageProperties
	QVTKImageWidget, 20
setAndDisplayMultipleImages, 20	
setImageProperties, 20	setImageStack
setXPicked, 20	ProbeCalibrationWidget, 15
setYPicked, 20	SetImageWidget
vtkImage, 22	QVTKImageWidgetCommand, 23
xPicked, 22	setMainWindow
xPosition, 22	ProbeCalibrationWidget, 15
yPicked, 22	SetPicker
yPosition, 22	QVTKImageWidgetCommand, 23
QVTKImageWidget.cpp, 27	setXPicked
- · · · · · · · · · · · · · · · · · · ·	QVTKImageWidget, 20
QVTKImageWidget.h, 28	setYPicked
ImageType, 28	
RGBImageType, 28	QVTKImageWidget, 20
RGBPixelType, 28	toutOn Tout A roo
QVTKImageWidgetCommand, 22	textOnTextArea
\sim QVTKImageWidgetCommand, 23	MainWindow, 12
Annotation, 23	transformations
Execute, 23	Calibration, 10
ImageWidget, 23	translations
New, 23	ProbeCalibrationWidget, 16
Picker, 23	
•	Ui, 7
QVTKImageWidgetCommand, 23	ui
QVTKImageWidgetCommand, 23	MainWindow, 13
SetAnnotation, 23	
SetImageWidget, 23	vtkImage
SetPicker, 23	QVTKImageWidget, 22
QVTKImageWidgetCommand.cpp, 29	arrimagornagoi, ==
QVTKImageWidgetCommand.h, 29	workWithStack
gvtkWidget	ProbeCalibrationWidget, 16
QVTKImageWidget, 21	Troboodiibration viagot, To
a v magomago, <u> </u>	xPicked
RGBImageType	QVTKImageWidget, 22
QVTKImageWidget.h, 28	
RGBPixelType	xPosition
QVTKImageWidget.h, 28	QVTKImageWidget, 22
	D: 1 1
renderer	yPicked
QVTKImageWidget, 21	QVTKImageWidget, 22
renwin	yPosition
QVTKImageWidget, 21	QVTKImageWidget, 22
rgbltklmage	
QVTKImageWidget, 22	
rotations	
ProbeCalibrationWidget, 16	